



THE UNIQUE DUAL INTER-LOCKING TABS SPIRAL AROUND THE NASOGASTRIC TUBE ENSURING SUPERIOR TUBE SECURITY.

The Dale NasoGastric Tube Holders are ideal for holding nasogastric feeding and aspiration tubes such as: Salem Sump® tubes, Levin stomach tubes, nasoenteric feeding tubes and naso-gastric intestinal tubes.

The holder is an improvement over tape and other types of adhesive backed holders, because it features a skin-friendly adhesive that adheres to the skin of the nose quickly and comfortably. The



holder can remain in place for up to 72 hours. The holder centers the tube exiting either from the left or right nostril, allowing the tube to "hang from" the nose in a more natural position.

- Easy to use.
- Securely holds all nasogastric tubes in place.
- Prevents nasal irritation, erosion or necrosis.
- Helps prevent inadvertent dislocation or removal of tube.
- Provides maximum tube security.
- Nose pad stretches and conforms to various nose shapes.
- Non-adherent tips allow holder to be removed easily without sticking to exam gloves.
- Not made with natural rubber latex.

| ORDERING INFORMATION | | |
|----------------------|--|------------|
| Product Number | Size | Quantity |
| 160 | NasoGastric Tube Holder Large/Adult | 50 per box |
| 161 | NasoGastric Tube Holder Small/Young Adult/Pediatric/Geriatric | 50 per box |

